


Antegrade Suture Passer



Antegrade Suture Passer



Middle-front jaw slot allows for easy loading

Notched, nitinol needle designed for reliable suture passing & security



Auto-capture window reliably traps and secures suture after passing

Large jaw opening of 15 mm in height and 18 mm in depth



Ergonomic handle designed for comfort and ease of use



www.maruho-medical.com

Suture Passer

Part #	Qty	Description	Outer Diameter	Working Length
26-20009	1	Suture Passer Instrument	4.5 mm	170 mm
26-30001	6	Disposable Needle	1.5 mm	240 mm

maruho
medical

200 Cobb Pkwy N
Building 200, Suite 210
Marietta GA 30062
888-404-3980
www.maruho-medical.com